



WHILL



Learn More Today at www.whill.inc

Model Ci2



WHILL Model Ci2

Adventure Awaits.


Where WHILL You Go?



Features and Functions



The battery can be removed and charged in a regular household outlet.



Effortlessly climb over obstacles up to 2" in height with two powerful motors and large front omni-wheels.



Adventure with WHILL - the Model Ci2 can make tight turns, allowing you to go places that a typical scooter can't take you.



Choose from 6 different arm color options and customize the Model Ci2 to match your personal style.



Dimensions: 17" (W) x 11.4" (L) x 9.5" (H)

The basket* is conveniently located under the seat, making it easy to get your items in and out.

*Optional feature



Remotely lock and unlock the Model Ci2 from a distance, with the smart key*.

*Optional feature

Accessories



Arm Covers
Choose from six colors.



Lap Belt
Firmly stabilizes the sitting position even on slopes.



Battery
The 2.7 kg lightweight battery can also be used as a spare battery.



Smart Key
Lock and unlock your Model Ci2 from a distance.



Battery Charger
A great compliment to a spare battery.



Dust Cover
Protects the Model Ci2 from dust and dirt.



Cane Holder
Canes and crutches can be fixed to the main body.



Armrest Pad
Easy to add cushion provides additional comfort and support.



Joystick
Operate joystick by grasping with your hand or holding it between fingers.

Specifications

Driving Range	11.2 mi / 18 km
Maximum Speed*	5 mph / 8 kmph
Obstacle Clearance	2"
Incline Capability	10°
Total Size (W x L x H)	21.8" x 38.8" x 29.3" - 37.2"
Weight	116 lbs
Charging Time	5 hr
Front Wheel	Omni wheel
Rear Wheel	Standard: Airless Tire
Arm	Lift-up type
Max Weight Capacity	300 lbs
Battery	Lithium ion (25.3V 10.6Ah)
Temperature Environment**	5 - 104°F

* The degree of deceleration may change when driving on steep slopes.
** If the battery temperature drops significantly while driving, the battery may go into cold mode and the deceleration rate may change.

WHILL Inc.

951 Mariners Island Blvd., Suite 300, San Mateo, CA 94404
(844) 699-4455 us.info@whill.inc www.whill.inc



The WHILL Model Ci2 is not considered to be a medical device and has not been submitted to the Food and Drug Administration for review or clearance.

